

ColorWorkDesk – LAB F TOP

FORMULATION AND CORRECTION SOFTWARE



“ColorWorkDesk LAB F Top”

The software in a single screen reports all the functions for the formulation, correction and storage of data.

The Navigator can manage standards and batches in a structure of folders and subfolders defined by the user to quickly and easily track information.

There are optimized formulation models for both covering and semi-covering products created by the ColorWorkDesk LAB M Top system and the formulation with container product line.

The active displays allow you to change the data directly where they are displayed and it is possible to save the work at each stage of the process.

The evaluation of the various formulas can be customized based on the customer's needs.

The automatic and manual corrections are very accurate, with the opportunity to ask the system to introduce a colorant foreign to the formula, to improve the DE in extremely critical situations.

The program is equipped with different types of graphs and prints for the evaluation of formulas, samples and colorants.

The ColorWorkDesk LAB F Top software does not include the program for creating and maintaining colorant files and product line.



FUNCTIONS AVAILABLE IN THE VERSION:

1. Color Scales: CIELab, CIELCh, XYZ, CIE_x, CIE_y, with absolutes and differences.
2. Quick storage of color charts in the database.
3. Management of gravimetric and volumetric units of measurement.
4. Creation of customized setups and masks for formulation and correction.
5. Color formulation in coverage, semi-coverage, transparency and with search & correct technology
6. Color formulation with product line
7. Automatic and manual SMART correction of the formula for the industry line.
8. Automatic and manual SMART correction of the formula for the product line.
9. Formula dosage with the creation of the output file to interface any dispenser.
10. Storage of formulas in progress and definitive.
11. Similar Search available directly in the formulation and correction module.
12. Formulation module in Queue with “background” work mode.
13. Export of results in txt and Excel format is guaranteed